



New

Rotosealer RS 120A

Advanced version Rotosealer

Offers practical means for controlling of all process parameters

- Improved temperature control
- Adjustable sealing pressure limit
- Microprocessor-controlled for reliable function
- All day ready for use
- TÜV tested
- Sealing technique according to DIN 58953-7
- Multi-language programming
- Audible alarm system upon incorrect values
- Simple control panel for easy use
- Exclusively Steriking®



Rotosealer RS 120A

A new model with advanced controlling and monitoring features

Offers practical means for controlling of all process parameters

- Provides improved safety and user convenience
- Automatic on/off function for reducing motor wear and for energy savings
- Built-in counter to calculate sealed packs
- Reliable operation for many years
- Simple control panel easy to use
- Adjustable in-feed guiding plate provides predefinition of the seal distance from the package edge



Technical data


Seal width	12 mm (2 x 5 mm)
Sealing speed	8,5 m/min
Sealing temperature	100-200 °C
Temperature tolerance	± 1%
Power consumption	180 W (max. 365W)
Power supply	220/240V 50/60 Hz 100/120 V 60Hz
Weight	12 kg
Dimensions	405 x 290 x 160 mm

Warranty

24 months from the date of purchase

Detailed operating manual and seal control sheet sample folder come with each unit.

Separate technical manual available on request.



We care that you pack safely!

▶ www.wipak.com

▶ www.steriking.fi

▶ e-Mail: steriking@wipak.com

Wipak Oy, P.O. Box 45, Wipaktie 2, FI-15561 Nastola, Finland, Tel. +358 (0)20 510 3333

