

Ideal for:

# INFU Packaging

Secondary films for FFS machines  
and multichamber primary webs



infusion liquids

STERILE 

## Overwrap films

TOP	Thickness	Structure	Open	Barrier	Description	Thick-ness	BOTTOM
PPP	100	PP	Tear	Standard	Double wound	100	PPP
OPP 2070	90	Laminated	Tear	Standard	Pouches	90	OPP 2070
OPP 20100	120	BoPP/PP				120	OPP 20100
OPP 2070	90	Laminated	Tear	Standard	FFS	100, 120,	PPP
OPP 20100	120	BoPP/PP				150, 200	
ESP 105	105	Laminated	Peel	Standard	FFS	150	PPP 150 Peel
		OPET/PP					

- For steam sterilization
- As overwrapping for packed liquids
- For materials with oxygen barrier properties, be in contact with our sales personnel



## Primary webs for multi-chamber applications

ITEM	Thickness	Structure	Description
PPP INFUPeel2	200	Polyolefin based film with DMF number	Excellent peel open properties in weak seal

Wipak reserves right to make changes at any time without notice.

WIPAK Oy  
 P.O. Box 45  
 Wipaktie 2  
 FI-15561 Nastola, Finland  
 Tel: +358-20510 311  
 Fax: +358-20510 3333  
 email: ind@wipak.com  
 www.wipak.com