



# WINDOWBAG

For surgical instrument sets





# WINDOWBAG - For surgical instrument sets

Windowbags are safe and reliable sterile barrier for instrument trays and kits. Window allows breathability and transparent film good visibility.

## Technical Details

- Multilayer PE film or reinforced PE/PA/PE film
- Printable coated Tyvek® window

## Applications

- For packaging of custom-assembled sets of instruments
- Other surgical items such as drapes, individual towels and tubing sets
- For EO and Irradiation sterilization

## Benefits

- Easy to use with clean peel
- Reliable barrier during transport, storage and handling
- Good clarity and transparency

DuPont™ and Tyvek® are trademarks of DuPont de Nemours, Inc.

For further Information please contact:  
medicaldevice@wiicare.com  
www.wiicare.com

Wiicare® is a registered trademark of Wihuri Packaging Oy.

 | **WICARE**

Wipak reserves right to make changes at any time without notice.