

 | **WIICARE**



Wiicare® FitTop™ Non-Forming Webs

ESP

Superior Film for Steam





ESP - Superior Film for Steam

ESP is specially developed for sterilization pouches giving clean peel but still strong seals when sealed with Wipak® steam papers. Tinted films helps visual check of the seal.

Technical Details

- Laminated BoPET/PP film
- Available in clear or blue
- 52 microns

Applications

- Sealable with steam papers
- Hospital pouches and reels
- For steam sterilization

Benefits

- High puncture resistance and excellent mechanical properties to keep package intact even with demanding applications
- Wide sealing window
- Good peel properties for clean and aseptic opening of peel packs

For further Information please contact:
medicaldevice@wiicare.com
www.wiicare.com

Wiicare® is a registered trademark of Wihuri Packaging Oy.

 **WICARE**

Wipak reserves right to make changes at any time without notice.