



WINDOWBAG G2

Patented sterile barrier system concept



WINDOWBAG - G2

- Excellent protection of sterile medical devices and/or surgical draping packs
- Superior breathability for Ethylene Oxide sterilization
- Most convenient aseptic opening and sterile presentation

Excellent protection

- Unique multilayer structure film is used with strong and flexible features
- The nonpeelable top of the bag prevents any seal damage during packing and transportation

Superior breathability

- Use of uncoated Tyvek® provides superior breathability and does not limit the air and gas penetration / removal by overall coating
- Reduced risk of broken seals during sterilization process due to the position of the non peelable function of Tyvek®

Most convenient aseptic opening and sterile presentation

- WindowBag G2 gives wide peelable film-to-film window for opening
- No contamination risk from loose particles of coating adhesive
- The strong peelable seal facilitate clean, fibre-free, peel opening allowing aseptic presentation of the packed product

DuPont™ and Tyvek® are trademarks of DuPont de Nemours, Inc.

For further Information please contact:
medicaldevice@wiicare.com
www.wiicare.com

Wiicare® is a registered trademark of Wihuri Packaging Oy.

WIICARE

Wipak reserves right to make changes at any time without notice.