

STERIKING®

WIICARE



Sterile Barrier Systems & Products

for Hospitals, Dental Practices & Laundries

See-Through Peel Pouches & Rolls

S - Flat Pouches

For packing individual instruments, smaller sets and other items. (STEAM, EO, FO)

Art.code	Size (mm)	Pieces/case
S17	50 x 200	3.600
S1	50 x 250	2.400
S24	75 x 150	2.400
S2	75 x 200	2.400
S23	75 x 230	1.800
S4	75 x 270	3.600
S18	75 x 300	4.200
S22	75 x 520	2.400
S25	100 x 150	2.400
S3	100 x 200	1.800
S5	100 x 270	3.000
S8	100 x 300	1.200
S19	100 x 350	1.200
S9	100 x 400	1.200
S12	100 x 570	1.800
S27	120 x 400	1.200
S15	150 x 200	1.200
S6	150 x 270	1.800
S20	150 x 300	1.800
S26	150 x 350	1.800
S13	150 x 400	1.200
S28	150 x 520	1.200
S34	160 x 600	1.200
S33	180 x 320	1.000
S7	205 x 270	600
S7A	205 x 340	1.000
S10	205 x 400	600
S10A	205 x 440	600
S35	205 x 600	1.200
S2070	205 x 700	600
S2244	220 x 440	600
S16	250 x 380	1.200
S11	250 x 500	600
S36	250 x 600	500
S37	250 x 700	600
S29	270 x 350	600
S30	270 x 440	600
S21	300 x 500	600
S14	300 x 570	600
S38	320 x 500	600
S31	420 x 500	500
S32	420 x 600	500



S-D - Flat Pouches

(STEAM, EO, FO)

Art.code	Size (mm)	Pcs Per Shipper Box	Pcs per dispenser
S1-D	50 x 250	1000	200
S13-D	150 x 400	1000	200
S16-D	250 x 380	600	200
S18-D	75 x 300	2000	200
S2-D	75 x 200	1000	200
S20-D	150 x 300	1000	200
S24-D	75 x 150	2000	200
S26-D	150 x 350	1000	200
S3-D	100 x 200	1000	200
S4-D	75 x 270	1000	200
S5-D	100 x 270	1000	200
S6-D	150 x 270	1000	200
S8-D	100 x 300	1000	200
S9-D	100 x 400	1000	200

SdV/SSdV - Pouches for robotic instruments

for packing instruments used in a robotic surgery, and for other long devices.



Art.code	Size (mm)	Pieces/case
SdV10	200 x 800	400
SdV20	250 x 900	400
SSdV10*	200 x 800	400
SSdV20*	250 x 875	400

*with adhesive closure



B - Pouches with Gusset

For more voluminous items. (STEAM, EO, FO)

Art.code	Size (mm)	Pieces/case
B35	100 x 50 x 300	500
B30	100 x 50 x 360	500
B31	150 x 50 x 400	500
B32	150 x 50 x 460	500
B36	200 x 55 x 400	250
B33	200 x 55 x 500	250
B37	250 x 65 x 480	250
B34	300 x 80 x 550	250
B38	420 x 100 x 570	250

See-Through Peel Pouches & Rolls

R - Flat Rolls (STEAM, EO, FO)



Art.code	Size (mm)	Pieces/case
R39	50 : 200 m	5
R40	75 : 200 m	4
R41	100 : 200 m	3
R125	125 : 200 m	2
R42	150 : 200 m	2
R43	200 : 200 m	1
R44	250 : 200 m	1
R45	300 : 200 m	1
R46	350 : 200 m	1
R47	400 : 200 m	1

Art.code	Size (mm)	Pieces/case
R39/100	50 : 100 m	8
R40/100	75 : 100 m	8
R41/100	100 : 100 m	4
R42/100	150 : 100 m	2
R43/100	200 : 100 m	2
R44/100	250 : 100 m	2
R45/100	300 : 100 m	1
R47/100	400 : 100 m	1
R49/100	500 : 100 m	1

Art.code	Size (mm)	Pieces/case
R39/31	50 : 31 m	10
R40/31	75 : 31 m	10
R41/31	100 : 31 m	10
R42/31	150 : 31 m	16
R43/31	200 : 31 m	6
R44/31	250 : 31 m	6
R45/31	300 : 31 m	6

RB - Rolls with Gusset (STEAM, EO, FO)



Art.code	Size (mm)	Pieces/case
RB50	75 x 25 : 100 m	4
RB51	100 x 50 : 100 m	3
RB52	150 x 50 : 100 m	2
RB53	200 x 55 : 100 m	1
RB54	250 x 65 : 100 m	1
RB55	300 x 80 : 100 m	1
RB56	350 x 80 : 100 m	1
RB57	400 x 80 : 100 m	1

RR - Reinforced Rolls (STEAM, EO, FO)



Art.code	Size (mm)	Pieces/case
RR41	100mm : 100m	2
RR43	200mm : 100m	1
RR45	300mm : 100m	1
RR47	400mm : 100m	1
RR49	500mm : 100m	1

SS - Standard Self Seal Pouches

with adhesive closure.



Art.code	Size (mm)	Pieces/case
SS10	60 x 250	1.000
SS1	90 x 200	1.000
SS2	90 x 250	1.000
SS3	90 x 570	1.000
SS4	130 x 270	1.000
SS4A	130 x 380	1.000
SS5	200 x 350	1.200
SS5A	190 x 330	1000
SS6	250 x 400	600
SS7	300 x 450	600

SST - Self Seal Plus Pouches



with an attached internal indicator strip.


Art.code	Size (mm)	Pieces/case
SST1	90 x 200	1.000
SST2	90 x 200	1.000
SST3A	100 x 640	1.000
SST4	130 x 270	1.000
SST4A	130 x 380	1.000
SST5A	190 x 330	1.000
SST6	250 x 400	600
SST7	300 x 450	600

SNW - TEX-Line Pouches



Peel pouches with nonwoven and film for packing heavy and bulky items. (STEAM, EO, FO)

Art.code	Size (mm)	Pieces/case
SNW2040	205 x 400	500
SNW2744	270 x 440	500
SNW3250	320 x 500	500
SNW4260	420 x 600	300
SNW5065	500 x 650	200
SNW5075	500 x 750	200
SNW5775	570 x 750	200



CB - Dust Cover Bags

for use after sterilization.

Art.code	Size (mm)	Pieces/case
CB1	150 x 270	1.000
CB2	205 x 300	1.000
CB3	250 x 400	1.000
CB4	300 x 500	500
CB5	400 x 550	500
CB6	420 x 700	500
CB7	500 x 800	300

SSDC - Self Seal Dust Cover Bags

with adhesive closure.

Art.code	Size (mm)	Pieces/case
SSDC10	150 x 270	500
SSDC20	200 x 350	500
SSDC30	250 x 400	500
SSDC40	300 x 550	200
SSDC60	420 x 700	200

LTS - Flat Pouches

CE MD

Made with DuPont™ Tyvek® for low temperature sterilization. Printed with hydrogen plasma process indicator.

Art.code	Size (mm)	Pieces/case
LTS7520	75 x 200	1.000
LTS1025	100 x 250	1.000
LTS1530	150 x 300	1.000
LTS1644	160 x 440	600
LTS1660	160 x 600	600
LTS2038	205 x 390	1.200
LTS2538	250 x 390	600
LTS2550	250 x 500	600

LTSS - Self Seal Pouches

CE MD

Made with DuPont™ Tyvek® for low temperature sterilization. Printed with hydrogen plasma process indicator.

Art.code	Size (mm)	Pieces/case
LTSS1	90 x 200	1.000
LTSS2	90 x 250	1.000
LTSS4	130 x 270	1.000
LTSS4A	130 x 380	1.000
LTSS5A	190 x 330	1.000
LTSS6	250 x 400	600
LTSS7	300 x 450	600

LTR - Flat Rolls

CE MD

Made with DuPont™ Tyvek® for low temperature sterilization. Printed with hydrogen plasma process indicator.

Art.code	Size (mm)	Pieces/case
LTR40	75 x 100	2
LTR41	100 x 100	2
LTR42	150 x 100	1
LTR43	200 x 100	1
LTR43A	225 x 100	1
LTR44	250 x 100	1
LTR45	300 x 100	1
LTR46	350 x 100	1
LTR47	400 x 100	1
LTR49	500 x 100	1

Wrapping Sheets

SMX - High Performance ProWraps

SMX-1 Blue Light

CE MD

For small and medium-size trays and sets.

Art.code	Size (mm)	Pieces/case
SMX1060	600 x 600	480
SMX1075	750 x 750	200
SMX1090	900 x 900	160
SMX1100	1.000 x 1.000	160
SMX1120	1.200 x 1.200	120
SMX1120/140	1.200 x 1.400	104
SMX1140	1.400 x 1.400	104

SMX-2 Blue Regular

CE MD

For medium and large-size trays and sets and for container inner wraps.

Art.code	Size (mm)	Pieces/case
SMX2060	600 x 600	480
SMX2075	750 x 750	200
SMX2090	900 x 900	120
SMX2100	1.000 x 1.000	120
SMX2120	1.200 x 1.200	120
SMX2120/140	1.200 x 1.400	104
SMX2140	1.400 x 1.400	104

SMX-2G Green Regular

CE MD

For medium and large-size trays and sets and for container inner wraps.

Art.code	Size (mm)	Pieces/case
SMX2G075	750 x 750	200
SMX2G090	900 x 900	120
SMX2G100	1.000 x 1.000	120
SMX2G120	1.200 x 1.200	120

SMX-3 Blue Heavy Duty

CE MD

Heavier structure for the most demanding applications.

Art.code	Size (mm)	Pieces/case
SMX3060	600 x 600	416
SMX3075	750 x 750	120
SMX3090	900 x 900	120
SMX3100	1.000 x 1.000	120
SMX3120	1.200 x 1.200	104

SMX-4 Blue Super Duty

CE MD

Especially high thickness for orthopedic sets and other heavy-duty uses.

Art.code	Size (mm)	Pieces/case
SMX4060	600 x 600	160
SMX4075	750 x 750	120
SMX4090	900 x 900	104
SMX4100	1.000 x 1.000	104
SMX4120	1.200 x 1.200	80

BW - BONDED WRAPS

CE MD

BW - Bonded Wrap 1 (37+47gsm)

Product Code	Size (mm)	Pieces/case
BW1060	600 x 600	288
BW1075	750 x 750	104
BW1090	900 x 900	104
BW1100	1.000 x 1.000	88
BW1120	1.200 x 1.200	72
BW1140	1.400 x 1.400	64

CE MD



BW - Bonded Wrap 2 (47+47gsm)

Product Code	Size (mm)	Pieces/case
BW2060	600 x 600	256
BW2075	750 x 750	96
BW2090	900 x 900	96
BW2100	1.000 x 1.000	88
BW2120	1.200 x 1.200	64
BW2140	1.400 x 1.400	56

CE MD

BW - Bonded Wrap 4 (75+47gsm)

Product Code	Size (mm)	Pieces/case
BW4060	600 x 600	208
BW4075	750 x 750	80
BW4090	900 x 900	80
BW4100	1.000 x 1.000	72
BW4120	1.200 x 1.200	56
BW4140	1.400 x 1.400	48



Wipak reserves right to make changes at any time without notice.

SPCG - Crepe Paper Green

CE MD

Special grade paper creped to give increased flexibility.

Art.code	Size (mm)	Pieces/case
SPC40G	400 x 400	500
SPC50G	500 x 500	500
SPC60G	600 x 600	500
SPC75G	750 x 750	250
SPC90G	900 x 900	250
SPC100G	1.000 x 1.000	250
SPC120G	1.200 x 1.200	100

NWG52 - Nonwoven 52

CE MD

For general use, for laundry and instrument sets.

Also available interleaved with SMX ProWrap grades.

Art.code	Size (mm)	Pieces/case
NWG52-60	600 x 600	500
NWG52-75	750 x 750	200
NWG52-90	900 x 900	200
NWG52-100	1.000 x 1.000	200
NWG52-120	1.200 x 1.200	100


Absorbent Trayliners

For use together with surgical instrument trays to avoid wet packs, protect the barrier and instruments.

TL-Standard	Size (mm)	Pieces/case
TL25	250 x 300	1.000
TL40	400 x 500	500
TL50	500 x 300	500
TL3040	300 x 400	500

TLWE-Trayliners with edge		
TLWE2525	250 x 250	500
TLWE4825	480 x 250	500
TLWE5030	500 x 300	500

TLNW-Trayliners nonwoven		
TLNW1525	150 x 250	3.200
TLNW3025	300 x 250	1.200
TLNW3040	300 x 400	1.200
TLNW3050	300 x 500	600



Steriking® is registered trademark of Wipak.

STERIKING COMBI SETS



STERIKING COMBI SET CW5 SMX1/NWL - COMBI WRAP LIGHT (37+52 gsm)		
Art.code	Size (mm)	Pieces/case
CW5075	750 x 750	240
CW5090	900 x 900	200
CW5100	1000 x 1000	144
CW5120	1200 x 1200	120
CW5120/140	1200 x 1400	120

STERIKING COMBI SET CW6 SMX2/NWL - COMBI WRAP REGULAR (47+52 gsm)		
Art.code	Size (mm)	Pieces/case
CW6075	750 x 750	240
CW6090	900 x 900	200
CW6100	1000 x 1000	144
CW6120	1200 x 1200	120

STERIKING COMBI SET CW7 SMX3/NWL - COMBI WRAP HEAVY DUTY (55+52 gsm)		
Art.code	Size (mm)	Pieces/case
CW7075	750 x 750	240
CW7090	900 x 900	200
CW7100	1000 x 1000	144
CW7120	1200 x 1200	120

STERIKING COMBI SET CW8 SMX4/NWL - COMBI WRAP SUPER DUTY (75+52 gsm)		
Art.code	Size (mm)	Pieces/case
CW8120	1200 x 1200	120

STERIKING COMBI SET CW10 SMX4/SMX2G (75+47 gsm)		
Art.code	Size (mm)	Pieces/case
CW10075	750 x 750	120
CW10090	900 x 900	104
CW10100	1000 x 1000	104
CW10120	1200 x 1200	80
CW10140	1400 x 1400	80



Indicator Strips

AC = steam, EO-gas, and

PL = low temperature plasma.

Art.code	Strips/dispenser	Pieces/case
AC6	250	
AC250/500	250	1.000
EO250	250	1.000
PL500		500

Daily Test Sheet

For Bowie & Dick Test with a cloth package.

Art.code	Pieces/case
DT50	50

Helix Challenge Test

For a charge control and as a steam penetration test for hollow instruments with small lumina.

Art.code	Strips/dispenser	Pieces/case
HX100	100	25
HX400	400	25

Tapes

Indicator Tapes, S = steam, H = hot-air indicator and L = latex free

Art.code	Width/length	Pieces/case
ITS19GL	19 mm x 50 m	48
ITS25GL	25 mm x 50 m	36
JTB19SG (B = blue)	19 mm x 50 m	48
JTB25SG (B = blue)	25 mm x 50 m	36
ITH19	19 mm x 50 m	48

Closing Tapes with no imprint.

Art.code	Width/length	Pieces/case
ITN19	19 mm x 50 m	48
ITN25	25 mm x 50 m	36
ITN50	50 mm x 50 m	24
JTB19N (B = blue)	19 mm x 50 m	48
JTB25N (B = blue)	25 mm x 50 m	36

Tape dispenser for one roll.

Art.code	Pieces/case
TDS25	1



Sealing Machines

Model: RS 120A Rotosealer

Model: RS 300 Rotosealer with printer

Model: IS 250 Impulse sealer with cutting device

Roll Holders in 2 different widths

• RH750, RH900



Seal Control Sheet

For operational qualification of the sealing process.

Art.code	Strips/dispenser	Pieces/case
SC250 (paper)	250	1.000
TSC200 (Tyvek®)	200	800

Multi Seal Test

For qualification and testing of the sealing process; visual check for seal quality and dye penetration test for seal integrity.

Art.code	Pieces/case
MS300	300

PB - Paper Bags

Heat-sealable lacquer top and steam indicator.

CE MD

Art.code	Size (mm)	Pieces/case
PB1	90 x 50 x 125	1.000
PB2	110 x 30 x 190	1.000
PB3	125 x 50 x 250	1.000
PB4	140 x 75 x 250	500
PB5	140 x 50 x 330	500
PB6	165 x 50 x 380	500
PB7	180 x 95 x 380	500
PB8	190 x 65 x 330	500
PB9	250 x 100 x 380	500
PB10	300 x 75 x 530	250
PB12	380 x 125 x 610	250

HR - Hot Air Rolls

CE MD

Special polyamide tubing for sterilization by dry heat.

Art.code	Width/length	Pieces/case
HR39	50 mm x 200 m	5
HR40	75 mm x 200 m	4
HR41	100 mm x 200 m	3
HR42	150 mm x 200 m	2
HR43	200 mm x 200 m	1
HR44	250 mm x 200 m	1

Registration Labels

For registering the sterilized items, enabling recall of all items if necessary. PP-based labels, resistant to disinfection, washing and sterilization.

Art.code		Pieces/case
TRLAB	Tray Label unprinted	1.000
TRLABPR	Tray Label printed	1.000

Paper label with adhesive backing, suitable as a bar code carrier.


Art.code		Pieces/case
TRTAG	Tray Tag	1.000

Non-toxic Markers

Waterproof, heat-resistant markers for sterilization packages.

Art.code	Color Pieces	Size (mm)	Pieces/dispenser
MAR1	Black	0.75	10
MAR1T	Black	1.00	10
MAR2	Red	0.75	10
MAR3	Blue	0.75	10





STERIKING®

Global Leader in Hospital Sterilization Packaging

Offering:

- Pouches & Rolls
- Wrapping Sheets
- Chemical Indicators & Tapes
- Seal Validation Products
- Sealers and Roll holders

We care that you pack safely!



Wipak Oy
www.steriking.com
e-mail: steriking@wipak.com

