



WINDOWBAGS

Clear Protection, Confident Sterility

WINDOWBAGS

Clear Protection, Confident Sterility

Windowbags are designed to ensure safe, sterile, and efficient packaging for medical procedure trays, kits, and surgical packs. They offer a combination of visibility, breathability, and ease of use and support various sterilization methods and aseptic presentation.

They save time compared to picking up individual items separately and reduce waste due to less packaging materials being used.

Type of bag	Film	Breathable web	Features
Windowbag STD	Recycle-ready PE or PE/PA/PE	Coated Tyvek® 1073B or 2FS™	Easy to use; the most common header bag solution*
Windowbag G2	Recycle-ready PE	Tyvek® 2FS™ uncoated providing superior breathability	Wide peelable film-to-film window for aseptic presentation
Windowbag G3	Recycle-ready PE	Tyvek® 2FS™ uncoated providing superior breathability	Conventional film-to-film peel opening of the top part of the bag

*12 sizes of PE 100 WB + Coated Tyvek® 1073B Windowbags available from stock

The traditional Windowbag is manufactured from polyethylene or reinforced nylon/polyethylene film, combined with printable coated Tyvek®. This is the most common solution in the market for header bags, offering reliable performance and a clean peel. Wiicare provides a PE-film solution combined with Tyvek® 1073B, available directly from stock in 12 different sizes.

Windowbag G2 and G3 represent the next generation of Windowbags. They offer superior breathability for ethylene oxide (EO) sterilization thanks to the use of uncoated Tyvek®, and are combined with recycle-ready polyethylene film to form a fully recyclable preformed sterile barrier system. Additionally, they eliminate the risk of contamination from loose coating adhesive particles.

- G2: Offers the most convenient aseptic opening and sterile presentation, featuring a large peelable film-to-film window at the bottom of the packaging.
- G3: Also features a film-to-film peel opening but conventionally located at the header part of the bag.

Applications

- For packaging of custom-assembled sets of instruments
- Other surgical items such as drapes, individual towels and tubing sets
- For EO and Irradiation sterilization

For further Information please contact:
medicaldevice@wiicare.com | www.wiicare.com

